

---

# Example script to load and plot ISIS neutron training course data on MnF2 using Horace

## Table of Contents

Get data into format that can be plotted and fitted with Horace routines .....	1
Take a cut that encompasses the entire data range to show (Q,E) map .....	2
Take a 1d cut and plot it: .....	3

Assumes that you have already installed Horace and SpinW!

===== Russell Ewings - 20/5/2020

## Get data into format that can be plotted and fitted with Horace routines

```
%Data file
spe_file='C:\Russell\Software\ExcitationPowderPublish\data_files
\MAR18301_Ei12.00meV.nxspe';

%If in nxspe format should not need "par" file, which specifies
detector
%positions, as already contained in the data file. Use blank string
here
par_file='';

%Name of your choice for Horace sqw file
sqw_mnf2='C:\Russell\Software\ExcitationPowderPublish\data_files
\MnF2.sqw';

%Tell Horace we are using a direct geometry spectrometer
emode=1;

%Set incident energy
efix=12.0;

%Make the sqw file that we then work with later
gen_sqw_powder_test (spe_file, par_file, sqw_mnf2, efix, emode);

%=====

-----
Calculating limits of data for 1 spe files...
Time to compute limits:
Elapsed time is 0.019997 seconds
```

*CPU time is 0.015625 seconds*

-----  
*Creating output sqw file:  
Time to read spe and detector data:  
Elapsed time is 0.015999 seconds  
CPU time is 0.015625 seconds*


*Calculating projections...  
Time to convert from spe to sqw data:  
Elapsed time is 0.0059967 seconds  
CPU time is 0.0625 seconds*

-----  
*Taking cut from sqw object...  
Have data from 83790 pixels - now processing data...  
Sorting pixel information for 83790 pixels*

## Take a cut that encompasses the entire data range to show (Q,E) map

Notice the non-standard (for Horace) method of **not** specifying a projection axis

```
mnf2_cut=cut_sqw(sqw_mnf2,0.03,0.12,'-nopix');  
%In the above, 0.03 specifies the bin width in |Q|, 0.12 specifies the  
  bin  
%width in energy transfer. '-nopix' means we don't bother retaining  
%detector pixel information (see Horace documentation for further  
  details)  
  
%Plot this:  
plot(mnf2_cut);  
  
%Use Horace commands to change colour scale and axes limits (see  
  Horace  
%manual for further plotting options)  
lz 0 80;%colour scale  
lx 0 4.5;%x-axis limits  
ly 0 10;%y-axis limits  
  
Taking cut from data in file C:\Russell\Software  
\ExcitationPowderPublish\data_files\MnF2.sqw...  
Step 1 of 1; Have read data for 83790 pixels -- now processing  
data... -----> retained 83790 pixels
```



## Take a 1d cut and plot it:

```
%Uses the same syntax as for a 2d slice


mnf2_cut2=cut_sqw(sqw_mnf2,[1.2,1.4],[1,0.12,9], '-nopix');
%Here we've specified to integrate the signal between 1.2<Q<1.4, and
make a
%cut along the energy axis from 1 to 9 in steps of 0.12meV

plot(mnf2_cut2)

%If all you wish to do is plot the data, this should be enough to get
you
%started.

Taking cut from data in file C:\Russell\Software
\ExcitationPowderPublish\data_files\MnF2.sqw...
Step 1 of 1; Have read data for 5671 pixels -- now processing
data... -----> retained 2103 pixels
```

C:\Russell\Software\ExcitationPowderPublish\data\_files\MnF2.sqw  
 $Q \leq 1.4$  in  $[Q, 0, 0]$ ,  $-\text{Inf} \leq \text{dummy} \leq \text{Inf}$  in  $[0, \text{dummy}, 0]$ ,  $-\text{Inf} \leq \text{dummy} \leq \text{Inf}$  in  $[0, 0, \text{dummy}]$   
E=0.94:0.12:9.1



Published with MATLAB® R2019b